



CS SURGICAL, INC.

662 WHITNEY DRIVE

SLIDELL, LA. 70461

PHONE: 985-781-8292

FAX: 985-781-8244


E-MAIL: cssurgicalinc@aol.com

www.cssurgical.com

Specialty Instruments
for the
UROLOGIST

"QUALITY INSTRUMENTS YOU CAN DEPEND ON"

UROLOGICAL RETRACTORS & SYSTEMS



A. CS1000-60 JORDAN-SIMPSON PERINEAL RETRACTOR SET

This set is a stable retractor system for anterior and posterior urethroplasty, and other perineal procedures.

THIS SET INCLUDES: 1 EA. CS1000-61
1 EA. CS1000-620
2 EA. CS1000-62
2 EA. CS1000-63
1 EA. CS1000-64
1 EA. CS1000-65
1 EA. CS1000-66
1 EA. CS1000-67

CS1000-61 JORDAN-SIMPSON SEGMENTED RING, SLOTTED, LARGE, 8 1/2" x 12" I.D.
CS1000-62 JORDAN-VERNICK 3 PRONG RETRACTOR, W/ROTATIONAL HANDLE
CS1000-63 JORDAN-VERNICK 4 PRONG RETRACTOR, W/ROTATIONAL HANDLE
CS1000-620 JORDAN-SIMPSON SINGLE PRONG RETRACOR
CS1000-64 JORDAN-SIMPSON RIGHT ANGLED RETRACTOR, 1"
CS1000-65 JORDAN-SIMPSON RIGHT ANGLED RETRACTOR, 1 1/2"
CS1000-66 JORDAN-RICHARDSON-SIMPSON RETRACTOR, 1" x 1 1/2"
CS1000-67 GOULET-SIMPSON RETRACTOR, SMALL, 1"

NEW ITEM

CS1000-66-FO JORDAN-RICH RETRACTOR 1" x 1 1/2" W/FIBER OPTIC LIGHT CARRIER, UP ANGLE
CS1000-66-FO-1 JORDAN-RICH RETRACTOR 1" x 1 1/2" W/FIBER OPTIC LIGHT CARRIER, DOWN ANGLE

B. CS1000-69 JORDAN-GELMAN-ALLIS FORCEPS W/SIMPSON HANDLE

This instrument is useful for grasping the dissected urethra, and can be attached to the ring using a tilt

CS1000-614 JORDAN-SIMPSON PEDIATRIC PERINEAL SET

This set is the smaller version of the Jordan-Simpson perineal set, that is a widely accepted retractor system for anterior and posterior urethroplasty procedures.

THIS SET INCLUDES: 1 EA CS1000-610 JORDAN-SIMPSON SEG RING, SMALL
6" x 9" I.D., SLOTTED
1 EA CS1000-611 JORDAN-RICH RETRACTOR 3/4" x 1"
1 EA CS1000-612 JORDAN-RICH RETRACTOR 3/8" x 3/4"
2 EA CS1000-613 JORDAN-SIMP 3 PRONG RETR, DELICATE

All the above items are available individually

The above retractor systems and individual retractors are compatible and adaptable to the BOOKWALTER™ and the SIMPSON RETRACTOR SYSTEMS.



CS7001 GELMAN VISUAL URETHRAL SOUND 20FR I.D

RETRACTORS FOR RETROPUBIC PROSTATECTOMY

The retractors listed below provide optimal exposure of the dorsal vein complex and posterior retraction of the bladder.



CS1000-40 YU-HOLTGREWE MALLEABLE PROSTATIC RETRACTOR, W/SIMPSON HANDLE, 2 1/2" x 7"



CS1000-43 CARROLL MALLEABLE PROSTATIC RETRACTOR, W/SIMPSON HANDLE

The above retractors will adapt to the BOOKWALTER™, and the SIMPSON RETRACTOR SYSTEMS.



CS1000-41 YU-HOLTGREWE MALLEABLE PROSTATIC RETRACTOR, W/BALFOUR HANDLE, 2 1/2" x 7"



CSU530 STAMEY DORSAL VEIN APICAL RETRACTOR

The retractor blades provide a wide space on both sides of dorsal vein complex for placing a controlling suture

UROLOGICAL RETRACTORS, RINGS, & SYSTEMS



CS1000-50 WALTHER-SIMPSON PERINEAL SET

This set provides optimal exposure for radical perineal prostatectomy.

This set includes: 1 EA. CS1000-51
2 EA. CS1000-52
1 EA. CS1000-53
1 EA. CS1000-54

**CS1000-51 WALTHER-SIMPSON PERINEAL RING, OVAL,
ANGLED 25 DEGREES**

CS1000-52 WALTHER-SIMPSON MALLEABLE RETR., 1"x 4", W/PRO EDGE


This blade is useful for lateral retraction

CS1000-53 WALTHER-SIMPSON MALLEABLE RETR., 1 1/2" x 4" W/PRO EDGE

This blade is used as the rectum or posterior retractor

CS1000-54 WALTHER-SIMPSON NOTCHED RETRACTOR, 1 1/2 x 4"

This blade can be used as an anterior retractor



CS1000-48 MALLEABLE RETRACTOR W/ LIP, 1/2"x 3"

CS1000-45 SINGLE HOOK MALLEABLE RETRACTOR, W/ BALL PROBE

CS1000-46 DOUBLE HOOK MALLEABLE RETRACTOR, 8MM, W/ BALL PROBE


The above retractor blades will adapt the BOOKWALTER™ and the SIMPSON RETRACTOR SYSTEMS.



CS1000-61 JORDAN-SIMPSON SEGMENTED RING, LARGE, 8 1/2" x 12" I.D.

CS1000-610 JORDAN-SIMPSON SEGMENTED RING, SMALL, 6" x 9" I.D.

This ring is excellent for perineal radical prostatectomy, anterior and posterior urethroplasty, uro-gyn procedures, and for vaginal surgery. The new small ring provides excellent exposure in young adults.



CS5200 BALFOUR RETRACTOR SET, 8 1/2" SPREAD

CS5200-A BALFOUR RETRACTOR SET, 10 1/2" SPREAD

CS5200-AA BALFOUR RETRACTOR SET, 14" SPREAD

EACH SET INCLUDES: TWO PAIRS OF LATERAL SNAP-ON WIRE BLADES 2 1/2" AND 3 1/2" DEEP, ONE CENTER BLADE 3 3/8" x 2 3/4"

CS5060 FOURTH BLADE ATTACHMENT, RIGHT HAND ASSEMBLY

CS5061 FOURTH BLADE ATTACHMENT, LEFT HAND ASSEMBLY

EACH ATTACHMENT INCLUDES THE CHOICE OF EITHER A 2" x 6 1/2", OR 2 1/2" x 6 1/2" MALLEABLE RETRACTOR.

EACH ATTACHMENT WILL FIT ANY STYLE BALFOUR RETRACTOR




"QUALITY INSTRUMENTS YOU CAN DEPEND ON"

INSTRUMENTS FOR PERINEAL PROSTATECTOMY



CS882-225 LOWSLEY TRACTOR, 11" STR
CS882-235 LOWSLEY TRACTOR, 13" STR
CS882-240 LOWSLEY TRACTOR, 13" CVD

For applying gentle traction to the prostate without trauma to the gland




- A. CS903-104 YOUNG PROSTATIC RETRACTOR 8 3/4", ANTERIOR, 1 1/2" x 2 1/2" BLADE
- B. CS903-112 YOUNG PROSTATIC RETRACTOR 8 3/4", BIFURACTED , 1 1/2" x 2 1/2" BLADE
- C. CS903-120 YOUNG PROSTATIC RETRACTOR 8 3/4", NOTCHED
- D. CS903-138 YOUNG PROSTATIC RETRACTOR 8 3/4", LATERAL, 7/16" x 2" BLADE
- E. CS903-146 YOUNG PROSTATIC RETRACTOR 8 3/4", IMPROVED



CS903-148 YOUNG PROSTATIC TRACTOR 8"


" QUALITY INSTRUMENTS YOU CAN DEPEND ON "

ASSORTED SCISSORS



- A. CS16-1590 SUPER-CUT JAMISON SCISSORS, 6 1/4", CVD
- B. CS5-244 THOREK SCISSORS, 7/14", FULL CURVE
- C. CS5-254 THOREK SCISSORS, 10", FULL CURVE
- D. CS5-243 JORGENSON SCISSORS, 9", STRONGLY CURVE BLADES
- E. CS25-1280 SATINSKY "S" SHAPED SCISSORS, 10", STRONG CURVE

ASSORTED NEEDLE HOLDERS




- A. CS01-099 MASSON-LUTHY N.H., BAYONET SHANK, 8"
- B. CS01-100 MASSON-LUTHY N.H., BAYONET SHANK, 10"
- C. CS01-106 BOZEMAN-WERTHEIM N.H., OFFSET SHANK, 8"
- CS01-107 BOZEMAN-WERTHEIM N.H., OFFSET SHANK, 9"
- CS01-108 BOZEMAN-WERTHEIM N.H., OFFSET SHANK, 10 1/4"
- CS01-109 BOZEMAN-WERTHEIM N.H., OFFSET SHANK, 12"
- D. CS01-101 BOZEMAN-FINOCHIETTO N.H., BAYONET SHANK, ANGLED JAW, 8"
- CS01-102 SAME AS ABOVE, BUT 10"
- E. CS01-114 STRATTE NEEDLE HOLDER 8"
- CS01-115 STRATTE NEEDLE HOLDER, 10"
- F. CS01-116 STRATTE NEEDLE HOLDER, REVERSE, 8"
- CS01-117 STRATTE NEEDLE HOLDER, REVERSE, 10"

"QUALITY INSTRUMENTS YOU CAN DEPEND ON"

RIGHT ANGLE CLAMPS



- CS7-255 GEMINI CLAMP, 5 1/2"**
- CS5-330 GEMINI CLAMP, 7"**
- CS7-257 GEMINI CLAMP, 8"**
- CS7-258 GEMINI CLAMP, 9"**
- CS7-259 GEMINI CLAMP, 11"**



- A. CS22-6570 MEEKER RIGHT ANGLE CLAMP, 7"**
- CS22-6571 MEEKER RIGHT ANGLE CLAMP, 11"**
- B. CS25-836 KANTROWITZ RIGHT ANGLE CLAMP, 7 1/2"**
- CS25-837 KANTROWITZ RIGHT ANGLE CLAMP, 9 1/2"**
- CS25-838 KANTROWITZ RIGHT ANGLE CLAMP, 11"**

- CS29-875 KALMAN-MIXTER CLAMP, 11", DELICATE**
- CS29-876 KALMAN-MIXTER CLAMP, 11", STD. JAW**




- CS2005 McDOUGAL CLAMP, 10 1/2", 2CM JAW LENGTH, LEFT**
- CS2006 McDOUGAL CLAMP, 10 1/2", 2CM JAW LENGTH, RIGHT**

The left clamp should be used if the surgeon is passing forehand to the left.

The right clamp should be used if the surgeon is passing backhand to the right.

" QUALITY INSTRUMENTS YOU CAN DEPEND ON "

CS SURGICAL, INC.



- A. CS22-6550 MIXTER "BABY" FORCEPS, 5", CVD
- B. CS22-6570 MIXTER FORCEPS, 5 1/4", DELICATE 90° JAW



- CS516-220A-TC BABCOCK FORCEPS, 8", ANGLED 45°
- CS516-218A-TC BABCOCK FORCEPS, 9", ANGLED 45°
- CS516-224A-TC BABCOCK FORCEPS, 12", ANGLED 45°



CS29-310TC BIRKLE DORSAL CLAMP/NEEDLE HOLDER, 11", W/T.C. INSERT

This instrument is designed to aid in the fixation and suturing of the D.V.C. during a retropubic radical prostatectomy.

" QUALITY INSTRUMENTS YOU CAN DEPEND ON "



CGS55-550 DeBAKEY TISSUE FORCEPS, 6 1/4", 1.5MM JAW
CGS55-555 DeBAKEY TISSUE FORCEPS, 8", 1.5MM JAW
CGS55-560 DeBAKEY TISSUE FORCEPS, 9 1/2", 1.5MM JAW
CGS55-565 DeBAKEY TISSUE FORCEPS, 12", 1.5MM JAW

CGS55-566 DeBAKEY TISSUE FORCEPS, 6 1/4", 2.0MM JAW
CGS55-567 DeBAKEY TISSUE FORCEPS, 8", 2.0MM JAW
CGS55-568 DeBAKEY TISSUE FORCEPS, 9 1/2", 2.0MM JAW
CGS55-569 DeBAKEY TISSUE FORCEPS, 12" 2.0MM JAW

CGS55-570 DeBAKEY TISSUE FORCEPS, 6 1/4", 2.5MM JAW
CGS55-575 DeBAKEY TISSUE FORCEPS, 8", 2.5MM JAW
CGS55-580 DeBAKEY TISSUE FORCEPS, 9 1/2", 2.5MM JAW
CGS55-590 DeBAKEY TISSUE FORCEPS, 12", 2.5MM JAW

CGS55-595 DeBAKEY TISSUE FORCEPS, 6 1/4", 3.5MM JAW
CGS55-600 DeBAKEY TISSUE FORCEPS, 8", 3.5MM JAW
CGS55-605 DeBAKEY TISSUE FORCEPS, 9 1/2", 3.5MM JAW
CGS55-610 DeBAKEY TISSUE FORCEPS, 12", 3.5MM JAW



CS34-4409 GIL-VERNET RENAL SINUS RETRACTOR, MALLEABLE SHAFT, 8MM WIDE
CS34-4411 GIL-VERNET RENAL SINUS RETRACTOR, MALLEABLE SHAFT, 11MM WIDE
CS34-4415 GIL-VERNET RENAL SINUS RETRACTOR, MALLEABLE SHAFT, 15MM WIDE
CS34-4417 GIL-VERNET RENAL SINUS RETRACTOR, MALLEABLE SHAFT, 20MM WIDE
CS34-4423 GIL-VERNET RENAL SINUS RETRACTOR, MALLEABLE SHAFT, 23MM WIDE



CSU540 DURRANI DORSAL VEIN COMPLEX LIGATION NEEDLE

The tip design and length provide for relatively simple and predictable passage of the ligature between the D.V.C. and the urethra.

SUTURE GUIDES



THESE INSTRUMENTS ASSIST IN REANASTOMOSING THE URETHRA AND BLADDER NECK AFTER RADICAL RETROPUBIC PROSTATECTOMY


- A. CSU500 ROTH URETHRAL SUTURE GUIDE, SET OF TWO, 24 & 28FR**
CSU501-24 ROTH URETHRAL SUTURE GUIDE, 24FR
CSU501-28 ROTH URETHRAL SUTURE GUIDE, 28FR

THIS INSTRUMENT FEATURES FOUR SUTURE GUIDE GROOVES ARRANGED AT 90° INTERVALS

- B. CSU515 ROTH GRIP-TIP™ URETRHAL SUTURE GUIDE, SET OF TWO, 24 & 28FR**
CSU515-24 ROTH GRIP-TIP™ URETRHAL SUTURE GUIDE, 24FR
CSU515-28 ROTH GRIP-TIP™ URETHTRAL SUTURE GUIDE, 28FR

THIS INSTRUMENTS FEATURES: THREE RETRACTABLE GRIPPING RODS WITH NON-SLIP ENDS, AND THREE SETS OF PAIRED SUTURE GUIDE GROOVES ARRANGED AT 120° NTERVALS

" QUALITY INSTRUMENTS YOU CAN DEPEND ON "



- CSU403-0 STAMEY NEEDLE, STRAIGHT TIP**
- CSU403-15 STAMEY NEEDLE, 15° ANGLED TIP**
- CSU403-30 STAMEY NEEDLE, 30° ANGLED TIP**


For surgically curable female urinary incontinence

NO SCAPEL VASECTOMY INSTRUMENTS




- CS382-200 NO SCAPEL VASECTOMY SET, INCLUDES 1 EA CS28-300 & 1 EA CS28-305**
- A. CS28-300 NO SCAPEL DISSECTING FORCEPS, 5" CVD, SHARP POINTS**
- B. CS28-305 NO SCAPEL LOOP FIXATION FORCEPS, 5"**

MEATUS AND INCONTINENCE CLAMPS



- A. CSU413 ZISPER MEATUS CLAMP**
- B. CSU411 BLACK MEATUS CLAMP**
- C. CSU412 HYAMS MEATUS CLAMP**



- A. CSU414-L CUNNINGHAM INCONTINENCE CLAMP, LARGE**
- B. CSU414-R CUNNINGHAM INCONTINENCE CLAMP, REGULAR**
- C. CSU414-J CUNNINGHAM INCONTINENCE CLAMP, JUVENILE**
- D. CSU415 BAUMRUCKER INCONTINENCE CLAMP**



CS SURGICAL, INC.

662 WHITNEY DRIVE
SLIDELL, LA. 70461

PHONE: 985-781-8292
FAX: 985-781-8244
cssurgicalinc@aol.com
www.cssurgical.com

ASK FOR THE FOLLOWING CATALOGS AVAILABLE FROM *CS SURGICAL*:

SIMPSON TABLE MOUNTED RETRACTOR SYSTEM



VALUE LINE
OPERATING ROOM QUALITY
SURGICAL INSTRUMENTS

“LIFETIME GURANTEE”
MADE OF GERMAN STEEL
2012 CATALOG & PRICE LIST

Specialty Instruments
for the
**COLON & RECTAL
SURGEON**

GREENWALD

TUNGSTEN CARBIDE
OPERATING ROOM QUALITY
SURGICAL INSTRUMENTS
2012 CATALOG & PRICE LIST

“QUALITY INSTRUMENTS YOU CAN DEPEND ON”